



SABAMEDICARE

Point-of-care diagnostics · Lebanon

PRODUCT CATALOGUE

GCare Lipid Panel Meter

Exclusive distributor in Lebanon

SabaMedicare SARL

For prices, demos and installation, contact us:

WhatsApp +961 71 490 154

info@sabamedicare.com

sabamedicare.com


Hazmieh, Baabda, Lebanon

GCare *Lipid*

Manage your cardiac risk with 6 parameter Test



LIPID



URIC ACID



GLUCOSE



User-Friendly, GCare Lipid



Portable printer



Self Optical Check



Bluetooth & Wi-Fi






Detachable strip guide
Easy to Clean

GCare *Lipid*

Manage your cardiac risk with 6 parameter Test

System performance

■ Accuracy



■ Precision

Classification	TC		TG		HDL	
	Level 1	Level 2	Level 1	Level 2	Level 1	Level 2
Average (mg/dL)	113	277	125	268	47	53
CV (%)	3.5%	2.4%	4.0%	3.2%	1.9%	3.1%

* Standard Reference : HITACHI 7180

GCare Lipid Analyzer

Product name	GCare Lipid
Memory	Store 1000 Cholesterol Data / Store 500 Glucose Data / Store 500 Uric Acid Data
Model	GCM-110, GCM-110W, GCM-110B
weight	120g
sample size	Cholesterol 40 μ l / Glucose 0.5 μ l / Uric Acid 15 μ l
Dimensions	130 x 65 x 30 mm3
Hematocrit	Cholesterol 30~50% / Glucose 20~60% / Uric Acid 30~50%
Measurement time	Cholesterol 150 seconds / Glucose 5 seconds / Uric Acid 150 seconds
Power	(AAA) 1.5V X 4
Measurement Range	TC 100~400mg/dL / TG 50~600 mg/dL / HDL 25~100mg/dL / Glucose 10~900 mg/dL / Uric Acid 3~20 mg/dL
Measurement temperature	Cholesterol 18~30 $^{\circ}$ C / Glucose 4~40 $^{\circ}$ C / Uric Acid 18~30 $^{\circ}$ C
Measurement Range Unit	mg/dL or mmol/L
Storage temperature	-10~60 $^{\circ}$ C (14 $^{\circ}$ F~140 $^{\circ}$ F)
Storage humidity	10~75% R.H.

GCM-110C

GCare Lipid Analyzer
Printer (Option)



How to use

Ordering Information

Product Name	Package	Order No.
GCare Lipid Analyzer	Unit	GCM110
GCare Lipid Profile Strip	10 Tests/Box	GLP110
GCare Lipid Profile C/S	2 Levels/Box	GLP110CS
GCare Glucose Strip	50 Tests/Box	GBP110
GCare Uric Acid Strip	10 Tests/Box	Coming Soon



GCMS

GC Medical Science Corp.
26, Mugeuk-ro 65beon-gil, Geumwang-eup, Eumseong-gun,
Chungcheongbuk-do, 27632 Republic of Korea

EC REP MT Promedt Consulting GmbH

Ernst-Heckel-Straße 7, 66386 St. Ingbert, Germany

Tel.: +82-31-267-8458
Fax: +82-31-260-9416
E-mail: gcms.overseas@greencross.com

0606-404-57-1

Made in Korea