



SABAMEDICARE

Point-of-care diagnostics · Lebanon

PRODUCT CATALOGUE

Hipro A1 Plus Automated Immunoassay System

Exclusive distributor in Lebanon

SabaMedicare SARL

For prices, demos and installation, contact us:

WhatsApp +961 71 490 154

info@sabamedicare.com

sabamedicare.com

Hazmieh, Baabda, Lebanon

Supplied in Lebanon by SabaMedicare SARL. Manufacturer specifications follow.

SIMPLE OPERATION



1. Collect sample



2. Insert into cuvette



3. Insert into channel and test automatically

ADVANTAGES

Flexible and efficient

One-step operation
Single test package

Accurate and reliable

Multi-methodologies
3-level Calibration System
Third-party control from Randox

Fully automatic testing

Intelligent temperature control
Closed test without pollution risk



DETAILS *High throughput with rapid tests*



High throughput

100 Tests/ H
Minimum test time: 2min
Emergency priority



Cuvettes discarded automatically



Test items identified automatically
Build-in standard curve

Technical Parameter
Model HP-AFS/1-PLUS
W×D×H 555mm×300mm×365mm
Weight 17kg
Power AC110~265V

TEST MENU

Infection

CRP
SAA
PCT
Hs-CRP
Hs-CRP+CRP

Renal Function

mAlb
Cys-C
NGAL
α1-MG
β2-MG
RBP
Cre

COVID-19

IgM / IgG Antibody
Antigen

Microelement

25-OH-VD
Ferritin

Diabetes

HbA1c (NGSP)

Cardiac

H- FABP
Myo
cTnI
HCY
BNP
NT-ProBNP

LP(a)

Lp-PLA2
CK-MB
D-Dimer
HDL-C
LDL-C
TG
TCHO

Immunoglobulin

IgA
IgE
IgM
IgG
Fertility
HCG

Tumor Marker

PSA
fPSA
CEA
AFP
CA199
CA125 ☆
CA153 ☆

Gastric Function

PGI
PGII
Hp

Rheumatism

ASO
RF
Anti-CCP

Liver Function

PA
TRF
ADA
C3
C4

☆ Under Developing

17

years
in IVD

3

R&D
centers

150

patents

100

countries

62

CE certificates



Technology and Innovation

- Reagent R&D Center in US;
- Instrument technology from Freiburg Germany;
- 150 international patents and 12 invention patents.

Production capacity

- 150000 Tests/day;
- Certified by TUV Rheinland with ISO 13485 and CE marks.

Facility

- World-class high capacity manufacturing facility (10000 sq.m.) with 1500 sq.m.current GMP plants and strict quality control system.
- A new industrial base was established in China Medical City, Jiangsu (35,000 sq.m.)

*A further step of
Immunoassay System*

Production base: No.3 Building,FangyiTech-park, 313zhujiangRd,Shijiazhuangcity, P.R.China

R&D Center :3938Trust Way, Hayward,CA94545,USA

Tel: + 86 311 83832628 Fax: + 86 311 83859087

Service Hotline: 400-0191-606

Email: IMD@hipro.us

Website: www.hipro.us



A1 PLUS

Automated Immunoassay System

Originate in Silicon Valley

High throughput

Full methodology

Automatic intelligent test

Hipro

Hipro Biotechnology Co.,Ltd

