



SABAMEDICARE

Point-of-care diagnostics · Lebanon

PRODUCT CATALOGUE

Hipro PalmF Fluorescence Analyzer

Exclusive distributor in Lebanon

SabaMedicare SARL

For prices, demos and installation, contact us:

WhatsApp +961 71 490 154

info@sabamedicare.com

sabamedicare.com

Hazmieh, Baabda, Lebanon

Supplied in Lebanon by SabaMedicare SARL. Manufacturer specifications follow.

Test Menu

Category	Abbr.	Assays
Diabetes	HbA1c	Glycosylated Hemoglobin
Thyroid	TSH	Thyroid Stimulating Hormone
	TT3	Total Triiodothyronine
	TT4	Total Thyroxine
	FT4	Free Thyroxine
	☆ FT3	Free Triiodothyronine
Fertility	FSH	Follicle-Stimulating Hormone
	PROG	Progesterone
	LH	Luteinizing Hormone
	AMH	Anti-Mullerian Hormone
	β-HCG	β-Human Chorionic Gonadotropin
	☆ E2	Estradiol
	☆ PRL	Prolactin
	☆ T	Testosterone
Infection	IL-6	Interleukin-6
	CRP	C- Reactive Protein
	HBP	Heparin Binding Protein
	PCT	Procalcitonin
Craniocerebral	S100-β	Central Nervous System Specific Protein
Cardiovascular	Hs-cTnI	High Sensitivity Cardiac Troponin I
	D-Dimer	D-Dimer
	ST2	Soluble Growth Stimulating Expression Gene 2
Microelement	25-OH-VD	25-hydroxy-Vitamin D
	VB12	Vitamin B12
	Folic Acid	Folic Acid
	☆ Fer	Ferritin

☆ Coming soon

R&D Center: 3938 Trust Way, Hayward, CA 94545, USA
Production Base: No.3 Building, Fangyi Science Park,
 No.365 Huai An Road, Shijiazhuang, Hebei, China.

Website: www.hipro.us
Email: hipro@hipro.us



PalmF

PalmF

Fluorescence Immunoassay Analyzer

Portable but Powerful
 Multi-scenario application

Rapid Detection
 Test in 155 seconds

Long-lasting Endurance
 Up to 584 consecutive tests

Hipro

Hipro Biotechnology Co., Ltd.



High Precision

30 times
of sensitivity



Wide Linearity

3 times of industry average

High Sensitivity

Unique optical path design

Dual Calibration


Optical path & curve calibration

Long-lasting Power

584 consecutive tests for
AAA batteries



Multi-scenario Application



Your powerful laboratory in hand

Smart System

Intelligent Interpretation

Avoid manual interpretation errors





Acousto-optic system

Avoid operation errors
Two-color indicator light
Smart power management



One-step Operation



Portability and Flexibility

Specimen types include serum, plasma and whole blood

Exquisite design, suitable for one-handed operation

Parameters

Name	PalmF Fluorescence Immunoassay Analyzer
Size	172 × 80 × 56mm
Weight	≤ 220g (exclude battery)
Power	≤1.5VA